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**Referencias:** 1. Vaishnavi A, *et al.* TRKing down an old oncogene in a new era of targeted therapy. *Cancer Discov.* 2015;5(1):25-34. 2. Lange AM, *et al.* Inhibiting TRK proteins in clinical cancer therapy. *Cancers (Basel).* 2018;10(4):pii:E105. 3. Amatu A, *et al.* NTRK gene fusions as novel targets of cancer therapy across multiple tumour types. *ESMO Open.* 2016;1(2):e000023. 4. Cocco E, *et al.* NTRK fusion-positive cancers and TRK inhibitor therapy. *Nat Rev Clin Oncol.* 2018;15(12):731-47. 5. Gatalica Z, *et al.* Molecular characterization of cancers with NTRK gene fusions. *Mod Pathol.* 2019;32(1):147-53. 6. Farago AF, *et al.* Clinicopathologic features of non-small-cell lung cancer harboring an NTRK gene fusion. *JCO Precis Oncol.* 2018;2018. 7. Sethi R, *et al.* Mammary analogue secretory carcinoma: update on a new diagnosis of salivary gland malignancy. *Laryngoscope.* 2014;124(1):188-195. 8. Okamura R, *et al.* Analysis of NTRK alterations in pan-cancer adult and pediatric malignancies: implications for NTRK-targeted therapeutics. *JCO Precis Oncol.* 2018;2018. 9. Khotskaya YB, *et al.* Targeting TRK family proteins in cancer. *Pharmacol Ther.* 2017;173:58-66. 10. de Lartigue J. TRK inhibitors advance rapidly in "tumoragnostic" paradigm. *Oncology Live.* 2017;18(15). <https://www.onclive.com/publications/Oncology-live/2017/vol-18-no-15/trk-inhibitors-advance-rapidly-in-tumoragnosticparadigm>. Accessed May, 2021. 11. Stransky N, *et al.* The landscape of kinase fusions in cancer. *Nat Commun.* 2014;5:4846. 12. Kumar-Sinha C, *et al.* Landscape of gene fusions in epithelial cancers: seq and ye shall find. *Genome Med.* 2015;7:129. 13. Marchiò C, *et al.* ESMO recommendations on the standard methods to detect NTRK fusions in daily practice and clinical research. *Ann Oncol.* 2019;30(9):1417-1427. 14. Garrido P, *et al.* Multidisciplinary consensus on optimising the detection of NTRK gene alterations in tumours. *Clin Transl Oncol.* 2021 Feb 23. doi: 10.1007/s12094-021-02558-0. 15. Hechtman JF, *et al.* Pan-Trk immunohistochemistry is an efficient and reliable screen for the detection of NTRK fusions. *Am J Surg Pathol.* 2017;41(11):1547-51. 16. Chen H, *et al.* Analysis of pre-analytic factors affecting the success of clinical next-generation sequencing of solid organ malignancies. *Cancers (Basel).* 2015;7(3):1699-1715. 17. Cleveland Clinic Laboratories. CC-SIGN® NTRK Plus Gene Fusion NGS Panel. <https://clevelandcliniclabs.com/2019/01/31/new-cc-sign-ntrk-gene-analysis/>. Accessed May, 2021. 18. Du X, *et al.* ALK-rearrangement in non-small-cell lung cancer (NSCLC). *Thorac Cancer.* 2018;9(4):423-430. 19. NeoGenomics. NTRK NGS Fusion Profile. <https://neogenomics.com/test-menu/ntrk-ngs-fusion-profile>. Accessed May, 2021. 20. Thermo Fisher Scientific. Oncomine NTRK Testing Flyer. 2018. Available for download at <https://www.oncomine.com/ngs-for-pathologists>. Accessed May, 2021. 21. Molecular MD. NGS Lung-Thyroid Fusion Cancer Panel. December 16, 2016. [www.molecularmd.com](http://www.molecularmd.com). Accessed May, 2021. 22. Penault-Llorca F, *et al.* Testing algorithm for identification of patients with TRK fusion cancer. *J Clin Pathol.* 2019;72(7):460-467. 23. Murphy DA, *et al.* Detecting gene rearrangements in patient populations through a 2-step diagnostic test comprised of rapid IHC enrichment followed by sensitive next-generation sequencing. *Appl Immunohistochem Mol Morphol.* 2017;25(7):513-523. 24. Su D, *et al.* High performance of targeted next generation sequencing on variance detection in clinical tumor specimens in comparison with current conventional methods. *J Exp Clin Cancer Res.* 2017;36(1):121. 25. Stack EC, *et al.* Multiplexed immunohistochemistry, imaging, and quantitation: a review, with an assessment of tyramide signal amplification, multispectral imaging and multiplex analysis. *Methods.* 2014;70(1):46-58. 26. Hung YP, *et al.* Evaluation of pan-TRK immunohistochemistry in infantile fibrosarcoma, lipofibromatosis-like neural tumour and histological mimics. *Histopathology.* 2018;73(4):634-644. 27. Rudzinski ER, *et al.* Pan-Trk immunohistochemistry identifies NTRK rearrangements in pediatric mesenchymal tumors. *Am J Surg Pathol.* 2018;42(7):927-935. 28. NeoGenomics. Pan-TRK. <https://neogenomics.com/test-menu/pan-trk>. Accessed May, 2021. 29. Kummur S, *et al.* TRK inhibition: a new tumor-agnostic treatment strategy. *Target Oncol.* 2018;13(5):545-556. 30. Luk PP, *et al.* Biomarkers for ALK and ROS1 in lung cancer: immunohistochemistry and fluorescent in situ hybridization. *Arch Pathol Lab Med.* 2018;142(8):922-928. 31. Empire Genomics. NTRK3 Break Apart FISH Probe. <https://www.empiregenomics.com/fish-probes/hd-probes/NTRK3+Break+Apart+FISH+Probe>. Accessed May, 2021. 32. Ali G, Bruno R, Savino M in lung cytology? *Arch Pathol Lab Med.* 2018;142(4):480-9. 33. Liu ML, *et al.* Rt-PCR gene expression profiling of RNA from paraffinembedded tissues prepared using a range of different fixatives and conditions. *Methods Mol Biol.* 2011;724:205-37. 34. Tatematsu T, *et al.* Investigation of neurotrophic tyrosine kinase receptor 1 fusions and neurotrophic tyrosine kinase receptor family expression in non-small-cell lung cancer and sensitivity to AZD7451 in vitro. *Mol Clin Oncol.* 2014;2(5):725-30. 35. Roche (via PR NewsWire). Roche launches first IVD pan-TRK immunohistochemistry assay. November 27, 2018. <https://www.prnewswire.com/news-releases/roche-launches-first-ivd-pan-trk-immunohistochemistry-assay-300755647.html>. Accessed May, 2021. 36. Roche Tissue Diagnostics. VENTANA pan-TRK (EPR17341) Assay: Package Insert 1017533EN Rev A. [https://productlibrary.ventana.com/ventana\\_portal/executeSearch.do](https://productlibrary.ventana.com/ventana_portal/executeSearch.do). Accessed May, 2021. 37. Solomon JP, Benayed R, Hechtman JF, Ladanyi M. Identifying patients with NTRK fusion cancer. *Ann Oncol.* 2019;30(Suppl8):viii16-viii22. 38. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Non-Small Cell Lung Cancer. V.2.2020. ©National Comprehensive Cancer Network, Inc 2019. All rights reserved. Published December 23, 2019. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 39. Roche (via PR NewsWire). Roche launches first IVD pan-TRK immunohistochemistry assay. Nov 27, 2018. <https://www.prnewswire.com/news-releases/roche-launches-first-ivd-pan-trk-immunohistochemistryassay-300755647.html>. Accessed May, 2021. 40. Roche Tissue Diagnostics. VENTANA pan-TRK (EPR17341) Assay: Package Insert 1017533EN Rev A. [https://productlibrary.ventana.com/ventana\\_portal/executeSearch.do](https://productlibrary.ventana.com/ventana_portal/executeSearch.do). Accessed May, 2021. 41. Fedchenko N, *et al.* Different approaches for interpretation and reporting of immunohistochemistry analysis results in the bone tissue - a review. *Diagn Pathol.* 2014;9:221. 42. Xu T, *et al.* Gene fusion in malignant glioma: an emerging target for next-generation personalized treatment. *Transl Oncol.* 2018;11(3):609-18. 43. Wang Z, Gerstein M, Snyder M. RNA-Seq: a revolutionary tool for transcriptomics. *Nat Rev Genet.* 2009;10(1):57-63. 44. Aguado C, *et al.* Fusion gene and splice variant analyses in liquid biopsies of lung cancer patients. *Transl Lung Cancer Res.* 2016;5(5):525-31. 45. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Rectal Cancer. V.1.2020. ©National Comprehensive Cancer Network, Inc 2019. All rights reserved. Published December 19, 2019. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 46. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Cutaneous Melanoma. V.1.2020. ©National Comprehensive Cancer Network, Inc 2019. All rights reserved. Published December 19, 2019. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 47. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Occult Primary (Cancer of Unknown Primary [CUP]). V.2.2020. ©National Comprehensive Cancer Network, Inc 2020. All rights reserved. Published January 31, 2020. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 48. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Colon Cancer. V.1.2020. ©National Comprehensive Cancer Network, Inc 2019. All rights reserved. Published December 19, 2019. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 49. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Head and Neck Cancers. V.3.2019. ©National Comprehensive Cancer Network, Inc 2019. All rights reserved. Published September 16, 2019. Accessed May, 2021. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org). 50. Yoshino T, *et al.* JSCO-ESMO-ASCO-JSMO-TOS: international expert consensus recommendations for tumour-agnostic treatments in patients with solid tumours with microsatellite instability or NTRK fusions. *Ann Oncol.* 2020;31(7):861-72.